

BF082 MIDDLEBROOK OADC GROWTH SUPPLEMENT

FORMULA

Ingredients:

(per vial sufficient for 450 ml medium)

Bovine albumin fraction V	2.50 g
Dextrose	1.00 g
Catalase	0.002 g
Oleic acid	0.025 g
Sodium chloride	0.0425 g
Distilled water	50.00 ml

Directions:

Warm up the refrigerated supplement to room temperature. Shake well to form uniform suspension. Aseptically add 1 vial to 450 ml of sterile, molten, cooled (45-50°C) Middlebrook 7H9 Agar Base (B735) / Middlebrook 7H10 Agar Base (B729) / Middlebrook 7H10 Agar Base, Special (B727)/ Middlebrook 7H11 Agar Base (B726) / Middlebrook 7H11 Agar Base w/o Malachite green (B578). Mix well and pour into sterile screw capped tubes or containers.

Precautions :

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: An antimicrobial supplement recommended for isolation and cultivation of Mycobacteria.

Storage: Between 2 to 8°C.

Packing: 1 vials