

TECHNICAL SHEET

BF083

NOVOBIOCIN SUPPLEMENT

FORMULA

Ingredients:

(per vial sufficient for 250/1000 ml medium)

Novobiocin	5.00 mg
------------	---------

Directions:

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water & aseptically add to 1000 ml of sterile, molten RS Medium Base (B388)/Modified Soyabean Bile Broth Base (B1490) for 250 ml. Mix well and pour into sterile Petri plates. It is also used in 250 ml of Tryptic Soya Broth, Modified (mTSB) (B1440).

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use : An antibiotic supplement recommended for the presumptive identification of *Aeromonas hydrophila* and detection of *E.coli* 0157:H7 from food and animal feeding stuffs.

Storage: Between 2 to 8°C.

Packing : 5 vials