

**BF088**

**PERFRINGENS SUPPLEMENT I**

**FORMULA**

**Ingredients :**

(Per vial, sufficient for 500 ml medium)

Sodium sulphadiazine            50mg

**Directions :**

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) Perfringens Agar Base O.P.S.P. (B671) along with sterile Rehydrated contents of 1 vial of Perfringens Supplement II (BF089). Mix well and pour into sterile petri plates.

**Precautions :**

1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :**An antibiotic supplement, recommended for the selective isolation of Clostridium perfringens .

**Storage:** Between 2 to 8°C.

**Packing :**5 vial