

**BF092**

**CETRINIX SUPPLEMENT**

**FORMULA**

**Ingredients :**

(per vial sufficient for 500 ml .medium)

Cetrimide	100.00 mg
Nalidixic acid	7.500 mg

**Directions :**

Rehydrate the contents of 1 vial aseptically with 2 ml sterile distilled water and add to 500 ml of sterile, molten Pseudomonas Agar Base (B695)/ Pseudomonas Agar Base (B1482). Mix well and pour Petri plates.

**Precautions :** 1. For Laboratory Use.  
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use :** An antibiotic supplement recommended for the selective isolation of Pseudomonas species

**Storage:** Between 2 to 8°C.

**Packing :** 5 vials