

**BF112                      CAMPYLOBACTER SELECTIVE SUPPLEMENT (KARMALI)**

**FORMULA**

**Ingredients:**

(Per vial, sufficient for 500 ml medium)

Sodium pyruvate	50.000 mg
Vancomycin	10.000 mg
Cefoperazone	16.000 mg
Cycloheximide	50.000 mg

**Directions:**

Rehydrate the contents of 1 vial aseptically with 5 ml of 95% ethanol. Mix well and aseptically add it to 500 ml. Sterile molten Karmali Campylobacter Agar Base (B1080) Mix well and pour into sterile Petriplates.

**Precautions :** 1. For Laboratory Use.  
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

**Use:** This supplement recommended for the selective isolation of thermo tolerant Campylobacter species.

**Storage:** Between 2 to 8°C.

**Packing:** 5 vials