

BF112 CAMPYLOBACTER SELECTIVE SUPPLEMENT (KARMALI)

FORMULA

Ingredients:

(Per vial, sufficient for 500 ml medium)

Sodium pyruvate	50.000 mg
Vancomycin	10.000 mg
Cefoperazone	16.000 mg
Amphotericin B	5.000 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 5 ml of 95% ethanol. Mix well and aseptically add it to 500 ml. Sterile molten Karmali Campylobacter Agar Base (B1080) Mix well and pour into sterile Petriplates.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: This supplement recommended for the selective isolation of thermo tolerant Campylobacter species.

Storage: Between 2 to 8°C.

Packing: 5 vials