

BF113 CAMPYLOBACTER SELECTIVE SUPPLEMENT W/HEMIN (KARMALI)

FORMULA

Ingredients :

(Per vial, sufficient for 500 ml medium)

Sodium pyruvate	50mg
Vancomycin	10 mg
Cefoperazone	16 mg
Amphotericin B	5 mg
Hemin	16 mg

Directions :

Rehydrate the contents of 1 vial aseptically with 5 ml. of 95% ethanol. Mix well and aseptically add it to 500 ml. Sterile molten Karmali Campylobacter Agar Base (B1080) Mix well and pour into sterile Petri plates.

Precautions : 1. For Laboratory Use.

2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: This supplement recommended for the selective isolation of thermo tolerant Campylobacter species.

Storage: Between 2 to 8°C.

Packing : 5 vial