

TECHNICAL SHEET

BF167 CC Supplement, Modified

FORMULA

Ingredients:

(Per vial sufficient for 500 ml medium)

| | |
|-----------------|-------|
| Amphotericin B | 5 mg |
| Chloramphenicol | 25 mg |

Directions:

Rehydrate the contents of 1 vial aseptically with 5ml 50% of ethanol. Mix well and aseptically add it to 500 ml sterile, molten, cooled (45-50°C) Sabouraud Dextrose Agar Base, Modified Emmons (B1110). Mix well and pour into sterile petri plates.

Precautions : 1. For Laboratory Use.
2. Follow proper, established laboratory procedures in handling and disposing of infectious materials.

Use: Recommended for the selective cultivation of pathogenic fungi.

Storage: Between 2 to 8°C.

Packing: 5 vials

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.