

BI391

SODIUM BISELENITE

Specifications:

CAS No.: 7782-82-3

Formula: HNaO3

Molecular Weight: 150.96

Appearance: White to cream free flowing powder

Solubility: Soluble in water

Precautions:

1. For Laboratory Use Only
2. Toxic / Dangerous for environment material. Handle carefully.

Storage : Below 30°C in cool, dry place.

Packing : 500 g

Refer disclaimer overleaf

Page 01 of 01

Disclaimer:

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related BIOMARKLABORATORIES publications.

The information contained in this publication is based on our in-house studies and market performance and is to the best of our knowledge true and accurate. BIOMARK LABORATORIES reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

BIOMARK Laboratories – INDIA

www.biomarklabs.com

TECHNICAL SHEET